# Nailit

Nail cleaning  
device

User Manual  
Manual del usuario  
Mode d'emploi  
Manuale d'uso  
Benutzerhandbuch

*TensCare™*

3

12

22

33

43

3

## Safety

The device should be used in strict compliance with the safety and operating instructions provided in this manual.

**Do NOT** use the device close to the womb during pregnancy.

Avoid direct eye contact from the laser radiation/blue light/infrared radiation to protect against the danger of excessive exposure to visible and infrared light.

**Caution:** Keep the device away from children. Children under the age of 18 years should only use the device under adult supervision.

**Do NOT** immerse your device in water or place it close to excessive heat such as a fireplace or radiant heater or sources of high humidity such as a nebulizer or kettle as this may cause it to cease to operate correctly.

**Do NOT** use this device if the clear plastic filter is broken or removed.

There are no user serviceable parts. No modification of this equipment is allowed.

4

**Do NOT** use the device for the following persons:

- Suffering from sunburn.
- Under medical care for diseases that involve photosensitivity.
- Receiving photosensitising medications.

If you experience any unexpected side effects such as itching, that occurs within the first 48 hours of using the device, medical advice should be sought prior to further exposure.

The device can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Children shall not play with the appliance and should not be left unsupervised.

5

## Operating Instructions

The device should be applied to each nail you would like to clean. Nailit can be used with both hands and feet.

1. Before using Nailit, be sure to thoroughly clean all nails which are to be used with the device.
2. Clamp the Nailit device onto the nail. Ensure the nail is facing the transparent, round window and therefore exposed to the blue light.
3. Press and hold the ON/OFF button for 3 seconds to switch on the device.
4. The ON/OFF button will flash green, indicating the session is underway. The nail is now exposed to the laser and blue light. The green light will continue to flash for the duration of the session (7 minutes).
5. At the end of the session, the green light will stop flashing and the device will switch off automatically.

This device has not demonstrated effectiveness in the treatment of fungal infections. It is cleared only for improvement in nail appearance.

## Frequency of Use

6

Each nail should be used with the device once a day. A visible improvement of the appearance of the nails is estimated to take 3-6 months when used daily. However, the effectiveness may vary between individuals. This device can be used at home.

## Charging

Nailit is powered by a rechargeable Li-ion battery. A USB-DC charging cable is included in the kit.

The battery should last 185 sessions. The battery has an estimated charging time of 2.5 – 3 hours from empty to fully charged. When the battery is running low, the LED indicator will flash red.

**To charge the battery:** Simply plug the charging cable in the device and charge with a standard USB charger. When charging the LED will flash green.

**Warning:** The battery in this device can only be replaced by qualified technicians.

7

The device cannot be operated while charging. When the battery is fully charged, the indicator will light green continuously.

**Warning:** Use only the charging cable supplied. Use of other cables could be hazardous and will negate the warranty.

## Cleaning

After every use, we recommend that you clean the inner rubber padding. Apply medical alcohol to a cloth and wipe clean. Do not share the device with other people without cleaning the rubber padding.

Cleaning and user maintenance must not be carried out by children without supervision.

## Troubleshooting

8

If your Nailit is not working properly, please check the following:

Problem	Possible causes	Solution
Battery light is red	Rechargeable battery has run out of power	Connect device to charger until LED continuously lights green. The device cannot be operated while charging
Device does not operate	Technical fault	Contact dealer/manufacturer

If problems persist, or to report unexpected operation or events, call TensCare or your local dealer for advice.

## Frequently Asked Questions

- 1. When will nail appearance improve?**  
Estimated time for visible improvement is 3 to 6 months. Effectiveness may vary between individuals.
- 2. Can the laser burn me?**  
No. Use of a Class 1M laser device should be pain-free with no side-effects or trauma.

## Technical Information

9

Laser type	1M (IEC 60825-1)
Laser wavelength	905nm
Blue light wavelength	470nm
Laser power	25W
Maximum pulse intensity	22W
Output plug	USB to DC (DC $\Phi$ 2.5x $\Phi$ 0.7x10 mm)
Battery	TW 902244 Li-ion battery 3.7V
Battery life	22 hours
Weight	50g
Dimensions	69.5 x 40.5 x 35.5 mm
Safety Classification	Internal power source. Designed for continuous use.

10

## Environmental Specifications:

Operating:	<b>Temperature range:</b> 5°C-40°C <b>Humidity:</b> 15-90% RH non-condensing
Storage:	<b>Atmospheric pressure:</b> 700hPa-1060hPa <b>Temperature range:</b> -10°C-60°C <b>Humidity:</b> 15% to 75% RH non-condensing <b>Atmospheric pressure:</b> 700hPa-1060hPa
Battery Specification	
Battery type	Lithium ion battery (Model: TW 902244)
Battery capacity	900mAh
Storage temperature	-10°C-60°C (≤ 1 month) -10°C-45°C (≤ 3 months) -10°C-28°C (≤ 1 year)
Storage humidity	45% - 75% RH
Charging Specification	<b>Current:</b> max 450mA <b>Voltage:</b> 5V <b>Power:</b> 2.5W

11

**Complies with:**

EN55015:2013 +A1:2015, EN IEC 61000-3-2:2019,  
EN61000-3-3:2013 +A1:2019, EN61547:2009,  
IEC 60825-1 :2014 (Third addition)

**Photobiological safety**

Exempt Group under EN62471:2006, does not pose any photobiological hazard.

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**TensCare Ltd.**